

DONECO

King Size Xpress Heat Heating Pad

DN0506



Instruction Manual

PRODUCT SPECIFICATIONS

Pad Size: 12" x 24"
 Power Supply: AC 120V, 60Hz, 180W
 Pad Material: Microplush (top); Polyester (bottom)
 Auto Shut-Off Time: 120 minutes

6 HEAT SETTINGS:

Warm-temperature range (~105°F)
 Low-temperature range (~112°F)
 Medium Low-temperature range (~119°F)
 Medium-temperature range (~126°F)
 Medium High-temperature range (~133°F)
 High-temperature range (~140°F)

PRODUCT FEATURES




- Fast heating:** Heats in seconds to provide fast pain relief to aching muscles.
- The extra large pad is sized at 12" x 24":** Perfect for relief on many different muscle groups.
- 6 Temperature Settings:** Digital LCD controller adjusts heat easily with 6 temperature settings to accommodate your comfort level.
- Microplush Material:** Thick microplush fibers offer maximum comfort for daily use.
- 2-Hour Auto Shut-Off:** Conserves energy and prevents burns caused by excessive heating.
- Machine Washable:** Helps maintain beauty and softness for years to come.

2

USING DN0506

Attach the Pad Connector to the Control Connector. Push both parts together securely to ensure a tight fit.



▲ NOTE: The Control Connector must be connected to the Pad Connector before plugging the unit into an electrical outlet.

- After connecting the Pad Connector the Control Connector, plug the unit a 120V AC electrical outlet.
- To turn the DN0506 on, press the power button. Use the Heat ▲ and Heat ▼ buttons to select the desired heat setting:

1=Warm 4=Medium
 2=Low 5=Medium High
 3=Medium Low 6=High

- To turn the DN0506 off, press the power button one more time.

▲ NOTE: If the unit does not turn on when you press the power button, make sure the Control Connector is securely plugged into the Pad Connector.

USING MOIST HEAT THERAPY

In some cases, you may find that moist heat therapy provides more comfort and tension relief. When using moist heat therapy, you should protect furniture and clothing from becoming damp.

3

- Dampen one side of the DN0506 with water by patting the surface with a wet sponge or cloth or use a spray bottle to apply a fine mist. Do not thoroughly soak the DN0506.
- Do not spray or dampen the connector.



DANGER: Do not immerse any part of the DN0506 in water or place under running water when plugged in. Doing so can damage the DN0506 and cause serious personal injury.



▲ NOTE: The Pad Connector and Control Connector must be completely dry before establishing an electrical connection.

AUTO-OFF FEATURE

To conserve energy and prevent skin burns caused by excessive use, the DN0506 includes an auto-off feature. The unit will turn off automatically after two hours of continuous use. You can turn this feature off by pressing the Auto Off button until the Stay On indicator light turns red. The unit will not turn off automatically after two hours of continuous use when the Stay On light is red.

▲ NOTE: After a power outage, the DN0506 automatically resets the auto-off feature.

WASHING INSTRUCTIONS

Do Not Dry Clean! Do Not Use Bleach! Do Not Use Wringer!

- Disconnect the DN0506 from the electrical outlet.
- Detach the Pad Connector from the Control Connector.
- Machine wash the DN0506 in cold water on the gentle cycle or hand wash. Line dry or tumble dry on high heat for five minutes.
- If using an automatic dryer, remove while damp and drape to finish drying.
- Do not iron the heating pad.
- Do not plug in the DN0506 until the Pad Connector is completely dry.

4

WARRANTY

The DN0506 comes with an industry-leading 5 Year Warranty that begins on the date of purchase.

The warranty applies to the heating pad and necessary parts and labor relating thereto. The warranty does not apply to damage resulting from failure to follow the operating instructions, accidents, abuse, alterations, or disassembly by unauthorized individuals.

ARE YOU 100% SATISFIED?

If you have product questions or concerns, don't hesitate to contact us directly at support@myfrieq.com

DONECO has an award-winning warranty, exchange, and customer service program that guarantees hassle-free solutions to any issue you might have within 24 hours!

5

IMPORTANT SAFETY INSTRUCTIONS DO NOT DESTROY

DANGER: TO REDUCE THE RISK OF BURNS, ELECTRIC SHOCK, FIRE, AND PERSONAL INJURY, THIS PRODUCT MUST BE USED IN ACCORDANCE WITH THE FOLLOWING INSTRUCTIONS:

- READ ALL INSTRUCTIONS.
- DO NOT USE WHILE SLEEPING.
- DO NOT USE ON AN INFANT OR ON AN ANIMAL.
- THIS PRODUCT IS NOT TO BE USED ON OR BY AN INVALID, A SLEEPING OR UNCONSCIOUS PERSON, A PERSON WITH DIABETES, OR A PERSON WITH POOR BLOOD CIRCULATION UNLESS CAREFULLY ATTENDED.
- DO NOT USE ON AREAS OF THE SKIN THAT ARE INSENSITIVE TO TOUCH OR HEAT.
- BURNS CAN OCCUR REGARDLESS OF CONTROL SETTING. CHECK SKIN UNDER PRODUCT FREQUENTLY.
- DO NOT USE IN AN OXYGEN-FILLED ATMOSPHERE.
- DO NOT USE PINS OR OTHER METALLIC OBJECTS TO FASTEN THIS PRODUCT INTO PLACE.
- DO NOT SIT ON OR AGAINST THE PRODUCT OR CRUSH IT. AVOID SHARP FOLDS IN THE FABRIC.
- DO NOT USE THE CORD AS A HANDLE.
- CHILDREN SHOULD NOT USE THIS UNIT WITHOUT ADULT SUPERVISION.
- PLACE THE UNIT ON TOP OF AND NOT UNDER THE PART OF THE BODY THAT REQUIRES HEAT THERAPY.

6

- ONLY USE THIS PRODUCT ON A 110-120V AC CIRCUIT.
- UNPLUG WHEN NOT IN USE.
- DO NOT USE WITH GENERATORS, POWER CONVERTERS, OR INVERTERS.
- DO NOT TAMPER WITH THIS PRODUCT IN ANY WAY. THERE ARE NO USER-SERVICEABLE PARTS. IF FOR ANY REASON THIS PRODUCT DOES NOT FUNCTION PROPERLY, CONTACT PURE ENRICHMENT AT SUPPORT@MYFRIEQ.COM
- DO NOT USE THIS PRODUCT WITH LINIMENT, SALVE OR OINTMENT PREPARATIONS AS SKIN BURNS CAN RESULT.
- NEVER PULL THIS PAD BY THE SUPPLY CORD.
- DO NOT BEND OR PINCH THE CORD.
- DO NOT USE THE PRODUCT IF THE CORD IS DAMAGED.
- NEVER LEAVE THE PRODUCT UNATTENDED, ESPECIALLY IF CHILDREN ARE PRESENT.
- LOOP CORD LOOSELY WHEN STORING. TIGHT WRAPPING CAN DAMAGE CORD AND INTERNAL PARTS.
- THIS PRODUCT HAS A POLARIZED PLUG (ONE BLADE IS WIDER THAN THE OTHER). TO REDUCE THE RISK OF ELECTRIC SHOCK, THIS PLUG IS INTENDED TO FIT IN A POLARIZED OUTLET ONLY ONE WAY. IF THE PLUG DOES NOT FIT FULLY IN THE OUTLET, REVERSE THE PLUG. IF IT STILL DOES NOT FIT, CONTACT A QUALIFIED ELECTRICIAN TO INSTALL THE PROPER OUTLET.

DO NOT ATTEMPT TO MODIFY THE PLUG OR DISABLE THIS SAFETY FEATURE IN ANY WAY.

- SAVE THESE INSTRUCTIONS.

7

DONECO

Tel: 1-909-595-1200 (US)
 (Mon.9:00-17:00 PST)